

WITec
focus innovations

TESCAN
PERFORMANCE IN NANOSPACE



Raman Imaging and Scanning Electron Microscopy

RISE Microscopy

Correlative microscopy on a new level –
Complementing ultra-structural SEM with
molecular Raman imaging

www.witec.de

www.tescan.com

RISE Microscopy

World's first fully-integrated Raman Imaging and Scanning Electron Microscope



Correlative Scanning Electron and confocal Raman Imaging for comprehensive sample analysis

A new dimension in imaging: ultra-structural SEM complemented with chemical compound information and molecular Raman imaging



RISE Microscopy is suitable for you...


...as a Raman newcomer you will benefit from the ease-of-use and intuitive measurement procedure

...as an experienced user you will appreciate the exceptional correlative microscope performance with the advantages of both techniques included in one instrument




Materials science, nanotechnology, polymer science, geosciences, life science, pharmaceutical industry etc. ...

...regardless of the field in which you are working, RISE Microscopy will empower you with its unique imaging capabilities.



a) Backscattered SEM image of a geological sample (diorite). **b)** SEM image overlaid with the Raman image. The different colors in the Raman image illustrate the various molecular compounds. Raman image: 100 µm x 100 µm, 300 x 300 pixels = 90,000 spectra, integration time: 34ms/spectrum. **c)** The corresponding color-coded Raman spectra display each molecular component of the sample. Sample courtesy of Dr. Christine Heim, Geowissenschaftliches Zentrum, University of Goettingen, Germany



a) SEM image of a graphene sample **b)** Color-coded confocal Raman image. The colors display the graphene layers and wrinkles. Image parameters: 20 µm x 20 µm, 150 x 150 pixels = 22,500 spectra, integration time: 0.05 s/spectrum **c)** SEM image overlaid with the confocal Raman image.

Scanning Electron Microscopy


The scanning electron microscope (SEM) is a powerful tool for achieving an outstanding level of magnification and resolution of nanometers. This gate to the nanoworld utilizes focused beams of electrons to reveal information about e.g. the morphology, topography or chemical composition of the sample. Furthermore, the performance of SEMs can be expanded by the combination of electron and focused ion beam source (FIB) to obtain access to further analytical features.

All TESCAN instruments are equipped with modern optics incorporating a unique Wide Field Optics™ design with a proprietary Intermediate Lens (IML) in addition to a conventional objective lens for a variety of working and displaying modes. Real time In-Flight Beam Tracing™ which ensures precise, direct and continuous control of the beam current and the spot size. Moreover, TESCAN was the first to offer true live stereoscopic imaging using advanced 3D Beam Technology which opens up the micro and nano-world for an amazing 3D experience and 3D navigation.


TESCAN instruments

Thermal emission systems, LaB₆ systems, field emission systems and combined FIB systems

- Ultra-fast scanning system
- In-Flight Beam Tracing™
- Modern Optics
- A Wide Range of Display Modes
- Automated Procedures and User-Friendly Software
- In-Beam Secondary and Backscattered Electron Detectors



Inelastic scattering of light from a molecule.



Example Raman spectrum of cellulose.

The Raman Principle

- Raman spectroscopy is a well-established and non-destructive method to analyze the chemical composition of a sample.
- A Raman spectrum shows the energy shift of the excitation light (laser) as a result of inelastic scattering by the molecules in a sample.
- Each molecule and chemical compound results in a specific Raman spectrum and can be easily identified by its unique Raman 'fingerprint'.

Additional sample information from the Raman spectrum:

- a. Peak intensity: Quantity/amount of a specific compound
- b. Peak shift: Identification of stress and strain states
- c. Peak width: Degree of crystallinity
- d. Polarization state: Crystal symmetry and orientation



SCANNING ELECTRON MICROSCOPY

RISE Microscopy – The Instrument

The RISE Microscope combines all features of a stand-alone SEM and a top-class confocal Raman imaging microscope within one instrument.

- Quick and convenient switching between Raman and SEM measurement
- Automated sample transfer from one measuring position to the other
- Integrated software interface for user-friendly measurement control
- Correlation of the measurement results and image overlay
- No compromise in SEM and Raman imaging capabilities



Confocal Raman Imaging


- The WITec Confocal Raman Microscopy and Imaging System combines Raman spectroscopy with confocal microscopy and enables confocal Raman imaging with the information of a complete Raman spectrum at every image pixel and a lateral resolution at the diffraction limit (~200 nm).
- Confocal Raman Imaging with unprecedented performance in speed, sensitivity, and resolution
- Outstanding depth resolution ideally suited for 3D image generation and depth profiles
- Ultrahigh-throughput lens-based spectroscopic system for highest sensitivity and best performance in spectral resolution
- Ultra-fast Raman imaging option with only 0.76 ms integration time per spectrum
- Non-destructive imaging technique: no staining or fixation of the sample required

06 In-Beam SE and in-Beam BSE as well as in-chamber SE and BSE detectors.



05 The sample remains inside the vacuum chamber during the complete measurements to ensure a convenient work flow with ease-of-use.






01 High-performance SEM with Schottky field-emission gun achieves ultra-high resolution of 1 nm at 15 kV and 1,4 nm at 1 kV.




02 Fully integrated confocal Raman microscope with excellent imaging capabilities and outstanding performance in speed, sensitivity, and resolution.




03 The GM chamber was developed to extend the features of TESCAN SEMs with even better analytical potential. Large number of ports enables all detectors and techniques to point to the common analytical working point.

04 Once positioned on the scan table, the sample is automatically transferred and re-positioned between the measuring procedures.

APPLICATIONS



a) Raman image of a perforated structure in Si using a 532nm excitation laser. The image shows the intensity of the Si Raman signal. The holes in the structure appear dark. 10 μm x 10 μm , 200 x 200 pixels = 40,000 spectra, 34 ms integration time per spectrum. **b)** The Cross Section on the right shows the intensity along the green line. The holes with a spacing of approximately 360 nm are clearly resolved. **c)** Overlay of the SEM and the Raman image.



a) Raman-SEM image overlay of a hamster brain tissue sample. In the color-coded Raman image the white brain matter is shown in Green and the gray brain matter in Red. Raman Image: 100 μm x 100 μm , 300 x 300 pixels = 90,000 spectra, 50 ms integration time per spectrum. **b)** The corresponding Raman spectra reveal the different spectral characteristics of the white and gray brain matter.



Raman-SEM image overlay of a LT GaAs sample. In the color-coded Raman image the gold substrate (yellow color) can be clearly distinguished from GaAs (in Red) correlated with the ultra-structure of the sample. Raman Image: 50 μm x 50 μm , 300 x 300 pixels = 90,000 spectra, 34 ms integration time per spectrum.

3D RAMAN IMAGING WITH RISE 3D




Structured GaN layer used in semiconductor research: **a)** SEM image **b)** Color-coded Raman image, violet: GaN, green/yellow: fluorescence along the edges of the structure, red: stress in material **c)** Color-coded RISE 3D Image: 3D view of the 3D structure. The material characteristics and their distribution within the sample can be clearly distinguished. Scan-volume 60 x 15 x 20 μm^3 ; 180 x 45 x 20 pixels = 162,000 Raman spectra.

Sample courtesy of Dr. Eberhard Richter, Materials Technology Department, Ferdinand Braun Institute, Berlin, Germany

COMPARISON BETWEEN RISE MICROSCOPY AND ENERGY-DISPERSIVE X-RAY SPECTROSCOPY (EDX)

RISE Microscopy and EDX analysis of a geological sample: **a)** Overlaid SEM-EDX image: Three different element groups can be distinguished (Si, O: orange; Si, Al, Fe, Ca: purple; Na: green). **b)** Raman spectral imaging of same sample area (22,500 spectra; integration time: 0.08 s/spectrum): Three spectral clusters can be differentiated (Quartz, Epidote, Albite). The inserts show spectral variations through different mineral orientations. **c)** Overlaid color-coded Raman-SEM image shows the distribution of the molecular compounds (Yellow/Orange: Quartz in different orientations; Blue/Purple: Epidote in different orientations, Green: Albite).



WITec GmbH

WITec is the leading German manufacturer of confocal and scanning-probe microscopes for state-of-the-art Raman, Atomic Force (AFM), and Scanning Near-Field Optical Microscopy (SNOM). WITec's headquarters is located in Ulm, Germany, where all WITec products are developed and produced. Branch offices in USA, Japan, Singapore, and Spain ensure a worldwide sales and support network. From the company's founding in 1997, WITec has been distinguished by its innovative product portfolio and a microscope design that enables combinations of the various imaging techniques within one system.

An example of the company's breakthrough development is the world's first integrated Raman-AFM microscope. To this day, WITec's confocal microscopes are unrivaled in sensitivity, resolution and imaging capabilities. Significant innovation awards document WITec's enduring success and innovative strength.

TESCAN Brno, s.r.o.

The leading provider of scientific instrumentation is well known for its innovation and openness to work with researchers and customizing applications to fit specific analytical needs. TESCAN's product range includes thermal emission systems, LaB₆ systems, field emission systems, FIB and Plasma FIB systems.

The TESCAN brand has, within 23 years of its existence, built a formidable reputation; 1600+ SEMs installed in over 60 countries are a testament to TESCAN's first-class quality and proven technology.

TESCAN Brno, s.r.o. was established as subsidiary of a multi-national company TESCAN ORSAY HOLDING after the merger (2013) of Czech company TESCAN, a leading global supplier of SEMs and focused ion beam workstations, and the French company ORSAY PHYSICS, a world leader in customized Focused Ion Beam and Electron Beam technology.



WITec Headquarters in Ulm, Germany



TESCAN Headquarters in Brno, Czech Republic



TESCAN is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au

WITec Headquarters

WITec GmbH
Lise-Meitner-Straße 6 . D-89081 Ulm . Germany
Phone +49 (0) 731 140700 . Fax +49 (0) 731 14070200
info@WITec.de . www.WITec.de

TESCAN Headquarters

TESCAN Brno, s.r.o.
Libusina tr. 1 . 623 00 Brno – Kohoutovice . Czech Republic
Phone +420 530 353 211 . Fax +420 530 353 415
info@tescan.cz . www.tescan.com