



Advanced XRF

# ZSX PrimusIII+

X-ray Fluorescence Spectrometer

X-ray  
Fluorescence  
Spectrometer




Cu Ni Co Zn Ga Ge As Se Br Kr  
 Zn Ag Cd In Sn Sb Te Xe  
 Hf Ta W Re Os Ir Pt Au Hg Yb Al Pr  
 Rf Db Sg Bh Hs Mt Uuo  
 La Ce Pr Nd P Sm Eu Gd Dy Ho Er  
 Th Pa U Pu Am Cm Bk Cf Es Fm

# Reliable Production Control


## Features

### Principle

Excellent resolution with Wavelength Dispersive XRF optics



### WDX vs. EDX: Spectrum Comparison




High resolution and high precision


High sensitivity for light elements

### Simple & Safe Sample Preparation

Ideal Optics for Powder Analysis



Reduced risk of pressed pellets breaking and damaging spectrometer



Falling particulate matter can contaminate beam path, increased risk of pellets breaking


## Applications

### Analysis of Iron Ore

#### Unique Compton Scattering Ratio Method with Theoretical Alpha Corrections

Using the Compton scattering ratio with theoretical alpha corrections significantly improves the accuracy of the calibration curve for analysis of Fe in iron ore. This method is also applicable for other nonferrous metal ore and concentrate analysis.

T. Fe in Ore Calibration



Theoretical Alphas with Compton Ratio for Iron Ores

No.	Type	1st Comp.	2nd Co.	Coefficient	Edit
1	Aluminum/iron	Al	Fe	1.0190E-002	Insert
2	Aluminum/iron	CaO	Fe	3.1620E-002	Insert
3	Aluminum/iron	FeO	Fe	3.1620E-002	Insert
4	Aluminum/iron	Al2O3	Fe	1.7720E-002	Add
5	Aluminum/iron	TiO2	Fe	8.5200E-002	Delete
6	Aluminum/iron	HgO	Fe	1.6600E-002	Insert
7	Aluminum/iron	P	Fe	5.8800E-002	Insert
8	Aluminum/iron	S	Fe	7.4070E-002	Insert
9	Aluminum/iron	Cl2O	Fe	9.8800E-002	Insert
10	Aluminum/iron	Cl2O	Fe	1.7720E-002	Insert
11	Aluminum/iron	Fe	Fe	2.8100E-002	Insert


### ASTM C114 Analysis of Powder Portland Cement

ZSX Primus III+ was tested to determine if the instrument meets ASTM C114 qualification using Portland cement standards. The results demonstrate the spectrometer's ability to meet ASTM C114 requirements for the analysis of Portland cements.

Repeat analysis results (sample: SRM1889a)

Comp.	Certified value	10x repeat measurement		Difference between duplicates		Difference between average of pellets and certified value	
		1st pellet	2nd pellet	C114 criteria	Result	C114 criteria	Result
CaO	65.34	65.407	65.421	0.2	0.01	0.3	0.1
SiO <sub>2</sub>	20.66	20.739	20.775	0.16	0.04	0.2	0.1
Al <sub>2</sub> O <sub>3</sub>	3.89	3.854	3.861	0.2	0.01	0.2	0.03
Fe <sub>2</sub> O <sub>3</sub>	1.937	1.912	1.91	0.1	0.002	0.1	0.03
SO <sub>3</sub>	2.69	2.68	2.69	0.1	0.01	0.1	0.05
MgO	0.814	0.864	0.865	0.16	0.001	0.2	0.05
K <sub>2</sub> O	0.605	0.61	0.612	0.03	0.002	0.05	0.01
TiO <sub>2</sub>	0.227	0.225	0.223	0.02	0.002	0.03	0.003
Na <sub>2</sub> O	0.195	0.201	0.199	0.03	0.002	0.05	0.005
P <sub>2</sub> O <sub>5</sub>	0.11	0.113	0.112	0.03	0.001	0.03	0.002
Mn <sub>2</sub> O <sub>3</sub>	0.2588	0.2586	0.2592	0.03	0.001	0.03	0.002
Cl	0.0019	0.003	0.002	0.03	0.001	0.03	0.001
ZnO	0.0048	0.0042	0.004	0.03	0.000	0.03	0.001


SiO<sub>2</sub> in Rock Calibration



### Fusion Bead Rock Analysis

Fusion bead method is an effective sample preparation technique to eliminate errors caused by heterogeneous effects. The low dilution fusion technique enables the analyses of both major and minor elements. Unique fusion bead correction software assures accurate analysis of rocks.

La in Rock Calibration



# Optimized for Powder and Metal Analyses

The tube above optics in ZSX Primus III+ minimize powder scatter in the vacuum chamber and facilitate fast sample preparation by eliminating the need for binders when performing powder sample analysis.

The evacuation and vacuum leak rates can be switched between slow and fast to optimize sample throughput for powder and metal samples.

Advanced XRF

## ZSX Primus III+ X-ray Fluorescence Spectrometer



### X-ray Tube

#### X-ray Tube

The standard 3kW, Rh target, end-window tube efficiently generates fluorescent X-rays from samples for both heavy elements (Rh-K lines) and light elements (Rh-L lines).

#### Primary Beam Filter

Eliminates scattering peaks from the X-ray tube target (Thomson and Compton scattering). Reduces background, improving peak:background ratio.

Filter	Function and analyzing elements
Ni 400	Eliminates the scattering peaks from Rh target; Rh, Ru, Cd, Ag...
Ni 40	Reduces background; Pb, As, Mo, Zr, Nb...
Al 125	Reduces background; Ti, Co, Fe, Ce, Pr, Nd...
Al 25	Reduces background and reduces scattering from the Rh-target (L lines); Cl, Cd (L), Ag (L) ...



### Analyzing Crystal

#### 10 position crystal changer




The three standard analyzing crystals in the Primus III+ are capable of analyzing elements from O. The 10-position crystal changer makes it easy to put together the optimum combination from a wide range of optional analyzing crystals, creating a custom solution tailored to meet the needs of your specific applications.

Element	4	5	6	7	8	9	11	12	13	14	15	16	17	19	-	23	24	-	92	Characteristics		
Crystal	Be	B	C	N	O	F	Na	Mg	Al	Si	P	S	Cl	K	-	V	Cr	-	U			
Standard																						
RX25																						
PET																						
LiF(200)																						
RX40																						
RX61																						
Ge																					High sensitivity	
Optional																						
LiF(220)																						High resolution
RX35																						
RX45																						High sensitivity
RX75																						
RX80																						

### Sample Stage

#### High Precision Sample Positioning

The high precision positioning of the sample ensures that the distance between the sample surface and X-ray tube is kept constant. This is important for applications that require high precision, such as the analysis of alloys.




#### Special Optics Reduce Errors Caused by Curved Sample Surfaces

Perform high-precision analysis using a unique optical configuration designed to minimize errors caused by non-flat surfaces in samples such as fused beads and pressed pellets.

### Detectors

The system includes a scintillation counter (SC) for heavy element analysis and gas flow proportional counter (F-PC) with very thin polymer window for light element analysis are equipped. Automatic center wire cleaning keeps the F-PC's performance in peak condition with minimum maintenance.


X-ray pulse height distribution



### Support System

#### Easy to Clean Sample Chamber

In the event that powder particles contaminate the sample chamber, cleanup is easy. Powders in the evacuation system can be caught in the powder trap.



#### Auto Sample Changer

Standard system uses one tray (12 samples). Possible to set up to 48 samples when using four trays.



# Intuitive and Easy Operation

## Routine Analysis

### Analysis Package

The spectrometer can be supplied with two types of analysis packages. The pre-calibration package (optional) calibrates the spectrometer prior to shipment. With the pre-calibration package, the spectrometer will be ready to perform analyses using the drift correction samples after on-site installation. The application package (optional) contains both standard and drift correction samples, as well as analytical parameters and an installation manual, making instrument setup quick and painless for non-experts.


### Pre-Calibration Package

FeCoNi Alloys

#### Analyzed Elements and Concentration Range

Element	Range mass %	Element	Range mass %
Al	0.006 - 3.07	Ni	0.006 - 99.6
Si	0.007 - 4.06	Cu	0.007 - 32.09
P	0.002 - 0.03	Se	- 0.19
S	0.001 - 0.03	Zr	0.002 - 0.058
Ti	0.002 - 0.03	Nb	0.001 - 5.19
V	0.001 - 0.03	Mo	0.001 - 15.45
Cr	0.002 - 0.03	Sn	0.0007 - 0.014
Mn	0.0015 - 0.03	Ta	0.002 - 0.75
Fe	0.02 - 0.03	W	0.007 - 17.98
Co	0.0014 - 49.4	—	—


#### Relationship between Analyzed and Certified Values



### Application Package

Sample Type	Analysis Mode	Sample Prep.	Analyzing Elements
Low Alloy Steel, Stainless Steel	Calibration Method		Si, Mn, P, S, Ni, Cr, Mo, etc.
Special Steel, Nickel Alloy	FP Method	Surface Polishing	Mn, Si, Cr, Ni, Co, Mo, W, Fe, etc.
Brass, Lead Free-Cutting Brass			Cu, Fe, pb, Sn, Zn
Clay			
Silica Sand			SiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub> , Fe <sub>2</sub> O <sub>3</sub> , TiO <sub>2</sub> , MnO <sub>2</sub> , etc.
High Alumina			
Magnesia	Calibration Method	Fusion Method	SiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub> , Fe <sub>2</sub> O <sub>3</sub> , CaO, MgO, etc.
Chrome • Magnesia			SiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub> , Fe <sub>2</sub> O <sub>3</sub> , MgO, Cr <sub>2</sub> O <sub>3</sub> , etc.
Zircon • Zirconia			SiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub> , Fe <sub>2</sub> O <sub>3</sub> , MgO, ZrO <sub>2</sub> , etc.
Polyethylene		Hot press	F, Na, Mg, Al, Si, P, S, Cl, Ca, Ti, Cr, Fe, Zn

#### X-R control chart



### Tools for Statistical Process Control

The ZSX software has a wide array of functionality including statistical process control required for routine operation in production environments. Type standardization using the bias correction function makes it possible to analyze a range of alloys using standard-based calibration. Drift correction can be pre-programmed to run at preset intervals using the periodic analysis function. Statistical process control via X-R and other control charts are available to facilitate production and analysis control. Trend line graphs displaying drift correction coefficients simplify analysis control.

## Intuitive Analysis Window

### EZ Analysis

All routine operations are integrated in a new "EZ Analysis" menu. The intuitive user interface makes it easy to run daily operations from a single window.

#### Just 3 Steps

- Step 1**  
Choose sample position
- Step 2**  
Choose application
- Step 3**  
Click "Start"




### Advanced EZ Scan

- EZ Scan analysis (standardless, option) with no-fuss settings
- Example set-up for powder samples:
- Accurate results with sample information



- 1 Component type and balance
- 2 Sample film used
- 3 Matching library
- 4 Sample size (for correct thickness)
- 5 Binder information



### SQX Analysis

The optional "SQX" program is the best available tool for analyses under the FP (fundamental parameter) method. The scan-based, standardless analysis software utilizes a unique fixed angle measurement feature and provides accurate analysis results using the power of the full FP method, including matching library functionality and automatic overlap correction.

#### SQX Analysis Results


Sample: CCRMP Syenite SY-3

Major components											unit: mass %
Component	SiO <sub>2</sub>	Al <sub>2</sub> O <sub>3</sub>	Fe <sub>2</sub> O <sub>3</sub>	MnO	MgO	CaO	Na <sub>2</sub> O	K <sub>2</sub> O	TiO <sub>2</sub>	P <sub>2</sub> O <sub>5</sub>	
SQX value	59.23	11.72	6.20	0.31	2.74	8.25	4.16	4.16	0.16	0.57	
Std. value	59.68	11.76	6.49	0.32	2.67	8.25	4.12	4.23	0.15	0.54	


Minor elements											unit: mass %
Element	As	Ba	Ce	Co	Cr	Cu	Er	Ga	Gd	La	
SQX value	0.0027	0.043	0.193	0.0027	0.0005	0.0026	0.0061	0.0026	0.0094	0.112	
Std. value	0.0019	0.045	0.223	0.0009	0.0011	0.0017	0.0068	0.0027	0.0105	0.134	
Element	Nb	Nd	Ni	Pb	Pr	Rb	S	Sm	Sr	Th	
SQX value	0.0163	0.063	0.004	0.0151	0.0162	0.0188	0.037	0.0091	0.0287	0.0999	
Std. value	0.0148	0.067	0.001	0.0133	0.0223	0.0206	0.051	0.0109	0.0302	0.1003	
Element	U	V	Y	Yb	Zn	Zr					
SQX value	0.061	0.0016	0.0655	0.0051	0.0236	0.027					
Std. value	0.065	0.005	0.0718	0.0062	0.0244	0.032					

Matching library : CCRMP Syenite SY-2

#### Results screen of SQX



## Dimensions



unit:mm

## Installation Requirements

Power	3Phase 200/208V 40A 50/60Hz Single phase 100-240V 50/60Hz (PC)
Earth Grounding	Grounding resistance less than 30Ω, Independent
Cooling water	Quality: tap water or equivalent Temperature: 30°C or less Flow: 5L/min or more
Room temperature	15-28°C (daily variation less than ±2°C within the range)
Humidity	75% RH or less
Vibration	2m/s <sup>2</sup> or less (Lower than human sensitivity level)
Counter gas	P-10 gas (Ar 90% Methane 10% mixture) Gas Pressure 0.15MPa

## Specification

		ZSX Primus III+
X-ray Source	X-ray tube	End window type Rh Target 3kW
	High Voltage Generator	High Frequency Inverter type Maximum rating: 3kW, 60kV-100mA
Optics	Maximum Sample Size	ø51mm X 30mm (H)
	Primary Beam Filter	4 Filters
	Diaphragm	4 positions (ø35, 30, 20, 10mm)
	Crystal	10 crystal changer Standard: LiF (200), PET, RX25
Counting	Detector	Heavy Element: SC Light Element: F-PC

- The numeric values of performance indicated in this brochure are based on the test results at Rigaku. Rigaku does not warrant that the identical values can always be obtained regardless of different operational environments.
- Company names and product names in this catalog are trademarks of the companies and/or registered trademarks.

Specifications subject to change without notice.



**Rigaku Corporation** Tokyo Branch

4-14-4, Sendagaya, Shibuya-ku, Tokyo 151-0051, Japan

Phone:+81-3-3479-0618 Fax:+81-3-3479-6112 e-mail:rinttyo@rigaku.co.jp

[www.Rigaku.com](http://www.Rigaku.com)

### Rigaku Corporation Head Office

3-9-12, Matsubara-cho, Akishima-shi, Tokyo 196-8666, Japan  
Phone:+81-42-545-8189 Fax:+81-42-544-9223 e-mail:rinttyo@rigaku.co.jp

### Rigaku Corporation Osaka Factory

14-8, Akaoji-cho, Takatsuki-shi, Osaka 569-1146, Japan  
Phone:+81-72-693-7990 Fax:+81-72-693-6746 e-mail:rinttyo@rigaku.co.jp

### Rigaku Americas Corporation

9009 New Trails Drive The Woodlands, Texas 77381-5209, USA  
Phone:+1-281-362-2300 Fax:+1-281-364-3628 e-mail:info@rigaku.com

### Rigaku Beijing Corporation

2601A, Tengda Plaza, No.168, Xizhimenwai Avenue, Haidian District, Beijing 100044, P.R.China  
Phone:+86-010-8857-5768 Fax:+86-010-8857-5748 e-mail:info@rigaku.com.cn

### Rigaku Europe SE

Gross-Berliner Damm 151, 12487 Berlin, Germany  
Phone:+49-30-6264035-0 Fax:+49-30-6264035-10 e-mail:rehq@rigaku.co.jp



**Rigaku is proudly represented in Australia and New Zealand by**

**AXT Pty. Ltd.**

**1/3 Vuko Pl., Warriewood  
NSW 2102 Australia**

**T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365**

**W. [www.axt.com.au](http://www.axt.com.au) E. [info@axt.com.au](mailto:info@axt.com.au)**