



THE NEW GENERATION
OF SCIENTIFIC TRAINING LABORATORIES
FOR NANOTECHNOLOGY EDUCATION



NANOEDUCATOR II

NANOEDUCATOR II

NANOEDUCATOR II is a scientific measurement training complex based on a scanning probe microscope (SPM); suitable for both: scientific research as well as nanotechnology education



FEATURES:

- Improved quality of scanning results and high efficiency due to the use of a new digital controller
- Low noise closed-loop scanner 100x100x10 μm
- AFM / STM atomic lattice resolution
- Easy adjustments
- Robust design
- Remote control via the Internet

SCIENTIFIC MEASUREMENT TRAINING COMPLEX NANOEDUCATOR II CONSISTS OF:



- 1** Scanning Probe Microscope:
Atomic Force Microscopy (AFM),
Scanning Tunneling Microscopy (STM),
nanolithography



- 2** Digital controller



- 3** Electrochemical etch station
for probe production



- 4** Handbook, manual, demo
samples and accessories

NANOEDUCATOR II

NANOEDUCATOR II is a complete easy-to-use system that helps teachers to educate the next generation of researchers in nanoscience by means of thorough hands-on training in all important nanotechnology areas.

It is designed to capture the students' interest in science and train future nanotechnologists using both AFM and STM techniques.

Robust and foolproof, NANOEDUCATOR II is providing an interdisciplinary education with a broad understanding of different fields of nanoscience.



NANOEDUCATOR II provides:

- Individual approach to education process
- Remote control from the teacher's computer
- Interactive learning

THE TRAINING LABORATORIES

The classroom equipped with NANOEDUCATOR II can be used to demonstrate, image and measure basic nanoscience principles in next-generation curricula for chemistry, physics and biology as well as hands-on experience in investigation of the nanoobjects (nanoparticles) and nanostructures, lithography and nanomanipulations.

THE CLASSES OF FIRST-GENERATION NANOEDUCATOR:



Eindhoven University of Technology,
The Netherlands



Krasnoyarsk State Pedagogical University
by V.P. Astafiev, Russia

5 STEPS TO BRILLIANT RESULTS:



1 Make a tungsten wire probe or select the probe from the set provided



2 Insert the probe into the probe holder



3 Insert the sample




4 Approach the sample with the probe

5


Scan



STM image of graphite surface.
Image size:
800x800 nm



AFM image of calibration grating.
Image size:
70x70 μm




AFM image of human erythrocytes.
Image size:
35x35 μm

SOFTWARE:

Easy-to-use software runs under the operating Windows XP®; and Mac® OS. It performs 2D and 3D image processing and nano-scaled manipulations. Also it has the ability to share data via iPhone™ and iPad™.

The managing program provides the following procedures:

- Measurement preparation
- Set and view of parameters of current configuration of the measurement head and controller
- Set and view of parameters of planned measurements
- Sample surface imaging
- Control of the relative position of the sample and tip
- Approach management
- Scanning and spectroscopy measurement management



Mac®, iPhone™ and iPad™ are trademarks or registered trademarks of Apple Inc.; Windows XP® is the registered trademark of Microsoft Corporation

APPLICATIONS:



Biology
(cells, viruses, bacteria, DNA)



Material science
(metals, semiconductors, dielectrics, polymers, photovoltaic cells)




Data storage devices
(hard drives, CDs)



Micro- and nano-structures
(gratings, self-organizing systems)

ATOMIC LEVEL RESOLUTION

NANOEDUCATOR II has incorporated special electronics which enable the switch from conventional imaging to atomic level resolution by a simple mouse click.




STM image of atomic grating on graphite.
Image size:
3.5x3.5 nm



SPECIFICATIONS

Probe	Tungsten wire probes which can easily be made by students by the electrochemical etch station Piezo-tube AFM probe Piezo-tube STM probe Sharp tip for STM
Measuring head	Head for piezo-tube probes
Scanner	Piezo-tube closed-loop scanner
Scanning	By sample, range 100×100×10 μm
Sample positioned	Manual, range 5×5 mm
Sample weight	Up to 40 g
Sample size	Diameter up to 25 mm, Thickness up to 10 mm
Approach system	By sample, motorized, range 15 mm
Optical control	Embedded USB camera



SPM MODES

Atomic Force Microscopy (AFM):

- Topography imaging
- Phase imaging
- Force imaging
- Force spectroscopy
- AFM lithography

Scanning Tunneling Microscopy (STM):

- Constant Current
- Constant Height
- V(Z) spectroscopy
- I(V) spectroscopy



Digital controller



NANOEDUCATOR II
Dimensions: 170×150×150 mm

NANOEDUCATOR II



NT-MDT is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au

www.ntmdt.com