



In Situ GAS &
HEATING solution
for TEM platforms

 **climate**



www.denssolutions.com




Informaticalaan 12, 2628 ZD Delft



www.facebook.com/DENSsolutions



www.linkedin.com/company/denssolutions



The Climate In Situ TEM Gas & Heating Solution enables dynamic studies of specimen transformations in gas and elevated temperature at sub-Angstrom level. Climate is the only environmental solution in the market that allows full correlation of the structural and chemical data including reaction product analysis due to the integration with dedicated DENSSolutions Gas analyzer. Convert your high vacuum TEM from a static imaging tool into a real-world research laboratory, this will enable you to speed up your development of new catalysts or other energy relevant materials and techniques.




Climate holder

Typical applications



Heterogeneous catalysis



Nanomaterial growth & synthesis




Corrosion of metals & alloys



Green energy materials

Selected publications




Time resolved in situ TEM images demonstrate the shape evolution of a Pd nanocrystals under 1 bar O_2 at 200 °C and the rebuilt Pd equilibrium structure by calculations. Scale bar: 5 nm.

Shape evolutions of a metal nanocrystal in environmental conditions

The authors demonstrate an atomic scale TEM observation of shape evolutions of Pd nanocrystals under oxygen and hydrogen environment at atmospheric pressure. Combined with multi-scale structure reconstruction model calculations, the reshaping mechanism is fully understood. These results give a direct insight into the behavioural response of nanoparticles to a 'real' reactive pressure environment, which is likely to improve the understanding of solid-gas reaction during catalytic applications

Zhang, Xun, et al. "In situ TEM studies of the shape evolution of Pd nanocrystals under oxygen and hydrogen environments at atmospheric pressure." *Chemical Communications* 53.99 (2017): 13213-13216.




In-situ selected area electron diffraction (SAED) images.

The reaction kinetics of nanoparticles

Acquiring the kinetics of gas–nanoparticle fast reactions under ambient pressure is a challenge owing to the lack of appropriate in-situ techniques. Now an approach has been developed that integrates time-resolved in-situ electron diffraction and an atmospheric gas cell system in TEM, allowing quantitative structural information to be obtained under ambient pressure with millisecond time resolution. The ultrafast oxidation kinetics of Ni nanoparticles in oxygen was vividly obtained. This study gives new insights into Ni oxidation and paves the way to study the fast reaction kinetics of nanoparticles using ultrafast in-situ techniques.

Yu, Jian, et al. "Fast gas-solid reaction kinetics of nanoparticles unveiled by millisecond in-situ electron diffraction at ambient pressure." *Angewandte Chemie* (2018).



HRTEM images extracted from video sequences recorded under a mixture of $CH_4:H_2O$ at pressures equal to 10-2 mbar and 1.2 mbar respectively in a-b) and c-f). Initial Ir + Ir_2O_3 composite nanoparticle at 290°C (a- b) is gradually further oxidized (c-d) before and almost complete homogenization as Ir_2O_3 at 460°C (e) and its final complete reduction into metallic Ir at 485°C (f).


Structural & morphological evolution of catalyst

The mechanism of interfacial synergistic catalysis for supported metal catalysts has long been explored and investigated in several important heterogeneous catalytic processes. The modulation of metal–support interactions imposes a substantial influence on activity and selectivity of the catalytic reaction, due to the geometric/electronic structure of interfacial sites. Although great efforts have been made to validate the key role of interfacial sites in water-gas shift (WGS) over metal catalysts supported on reducible oxides, direct evidence at the atomic level is lacking and the mechanism of interfacial synergistic catalysis is still ambiguous. This work provides a fundamental understanding of interfacial synergistic catalysis toward the WGS reaction, which is constructive for the rational design and fabrication of high activity heterogeneous catalysts.

Xu, Ming, et al. "Insights into Interfacial Synergistic Catalysis over $Ni@TiO_{2-x}$ Catalyst toward Water-Gas Shift Reaction." *Journal of the American Chemical Society* (2018).

Why Climate?

1



Experimental preparation made easy

1. Easy sample preparation

Experienced Application Engineers will provide you with tips and tricks to avoid carbon contaminations from the beam for drop casted particles and a detailed workflow for FIB lamella preparation

2. Predefine the experimental conditions

Calorimetry and Mass spectrometry data produced by Climate system before TEM studies allows you to more easily define the optimum conditions for your experiment

2



State of art environmental control

1. Dynamic mixing

Using a specially designed and patented mixing valve instead of a premixing tank allows you to change the gas composition on the fly and to vary relative concentrations with dilutions of 0.1 %

2. Fast switching

Defined gas channel and minimal gas volume inside the Nano-Reactor enable to change the gas environment within seconds

3. Independent control

Gas composition, pressure and flow rate can be selected independently with the widest experimental range


4. Clean experiments

The modular design concept enables easy on-site exchange of all critical components without gluing or welding

5. Temperature accuracy and stability in gas environment

4-point probe heating provides the most accurate temperature with 0.01 °C stability, even during gas flow

3



High impact results

1. Data cross-correlation

Crystal structure, morphology, chemistry, thermodynamics and kinetics in one experiment

2. Product analysis

The optimized gas analyzer and defined gas channel of the Nano-Reactor guarantees an accurate analysis of the environment at any time with a sensitivity of 5 ppm

3. High Stability

Atomic resolution in TEM and STEM is routinely achievable in static and flow modes

4. Optimized analytical capabilities

The optimized design enables EELS and large solid angle EDS collection

Software for accurate environmental control

Climate Software

Advanced control and data analysis

Guided system preparation

- A step by step procedure for the preparation of the system with multiple safety features included in the software, ensures a peace of mind

Stimulus control

- The 'Direct control mode' allows each user to set all relevant parameters for the experiment such as pressure, gas composition and the flow rate in the Nano-Reactor
- The 'Flowsheet' mode allows the advanced user to set every component of the GSS independently

Automate your experiments

- Design your experimental workflow and obtain reproducible results with an embedded profile builder

Safe operation

- The gas flow in the Nano-Reactor is monitored and includes an automatic safety check to stop the flow if required
- Flammability limit check
- Automatic data logging during the experiment



Camera integration

True experimental integration

Easy & Fast experiments


- One SW to control in situ parameters & imaging

No mistakes or misinterpretation

- Synchronized in situ data with images at any frame rate

Save your time with your analysis

- Analyze images as a function of in situ parameters




System specifications

| JEOL | Climate Air | Climate G | Climate G+ |
|----------------------------------|--------------------------------------|--------------------------------------|--------------------------------------|
| Heating control | Four point probe resistive feedback | | |
| Temperature range | RT to 1000 °C | RT to 1000 °C | RT to 1000 °C |
| Temperature Stability | 0.01 °C | 0.01 °C | 0.01 °C |
| Pressure range | Ambient | 0 - 1000 mbar | 0 - 1000 mbar |
| Polepiece compatibility | All | All | All |
| Alpha tilt range | URP, FHP limited HRP, WGP ±20 deg | URP, FHP limited HRP, WGP ±20 deg | URP, FHP limited HRP, WGP ±20 deg |
| Modular design | Replaceable tubing and holder tip | | |
| Resolution | ≤100 pm* | ≤100 pm* | ≤100 pm* |
| Drift rate at in situ conditions | < 0.5 nm/min | < 0.5 nm/min | < 0.5 nm/min |
| Gas Mixing | N/A | Discrete | Continuous |
| Gas switching | N/A | < 60 s | < 1 s |
| Gas flow rate | Static | 0 – 3 mln/min | 0 – 3 mln/min |
| Micro-Calorimetry | V | V | V |
| Mass Spectrometer | V | V | V |

| Thermo Fisher Scientific | Climate Air | Climate G | Climate G+ |
|----------------------------------|---|---------------|---------------|
| Heating control | Four point probe resistive feedback | | |
| Temperature range | RT to 1000 °C | RT to 1000 °C | RT to 1000 °C |
| Temperature Stability | 0.01 °C | 0.01 °C | 0.01 °C |
| Pressure range | Ambient | 0 - 1000 mbar | 0 - 1000 mbar |
| Polepiece compatibility | BioTWIN, S-TWIN, TWIN, X-TWIN, C-TWIN | | |
| Alpha tilt range | BioTWIN, S-TWIN, TWIN: ±35 deg X-TWIN, C-TWIN: ±25 deg | | |
| Modular design | Replaceable tubing and holder tip | | |
| Resolution | ≤100 pm* | ≤100 pm* | ≤100 pm* |
| Drift rate at in situ conditions | < 0.5 nm/min | < 0.5 nm/min | < 0.5 nm/min |
| Gas Mixing | N/A | Discrete | Continuous |
| Gas switching | N/A | < 60 s | < 1 s |
| Gas flow rate | Static | 0 – 3 mln/min | 0 – 3 mln/min |
| Micro-Calorimetry | V | V | V |
| Mass Spectrometer | V | V | V |

*Listed specifications are dependent on microscope configuration

Complete 'plug & play' package



1. Climate Gas & Heating TEM specimen holder
2. Nano-Chips starter pack
3. Heating Control Unit
4. Laptop with Climate and Digiheater software
- 5a. GSS light (Climate G); 5b. GSS (Climate G+)
6. Optional: Gas Analyzer
7. Supporting tools



Service and Support:

| | |
|-----------------------------|---|
| Product Warranty | 24 months* with a possibility of extension |
| Product Liability Insurance | Every system includes product liability insurance to cover unlikely damages to the TEM or operators |
| Regulatory compliance | CE, RoHS, FCC |
| Mechanical compatibility | Approved by TEM manufacturers |
| Radiation safety | According to TEM manufacturers compliance regulations |
| Service | Dedicated Field Service Engineer |

*Depending on configuration








INNOVATIONS THAT MATTER

DENSsolutions

Informaticalaan 12
2628ZD Delft
The Netherlands
Phone: +31 (0) 153 030 214
Email: info@DENSsolutions.com
www.DENSsolutions.com

To find your local distributor, please visit
www.DENSsolutions.com/distributors



-  **wildfire** Heating
-  **lightning** Heating + Biasing
-  **climate** Gas + Heating
-  **stream** Liquid + Biasing or Heating



DENS Solutions is distributed in Australia
and New Zealand by AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au

Copyright © 2018 DENSsolutions B.V. All rights reserved.

Unless otherwise indicated, all materials on these pages are copyrighted by DENSsolutions B.V. All rights reserved. No part of these pages, either text or image may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, electronic, mechanical or otherwise, for reasons other than personal use, is strictly prohibited without prior written permission.