

HUNK



UNCHAINED  
LABS

## Make stability routine

What if you knew up front how a small tweak in formulation changed the stability of your biologic, or if it'll aggregate down the road? The only way to get your hands on that info is to know your  $\Delta G$ , and the HUNK makes measuring it totally routine. So you'll know how changes in protein constructs, concentrations, and formulations affect stability now, not later.

$\Delta G$  measure stability

$\Delta G_{\text{trend}}$  predict aggregation

$K_D$  small molecule affinity



## Why $\Delta G$ ?

Biologics are in equilibrium between their native and denatured states and  $\Delta G$  is the only way to know how much of each you've got at room temp. A high  $\Delta G$  means your protein's less denatured, which is a good thing. A low  $\Delta G$  means it's less stable and more denatured – a not so great thing. So figuring out if your biologic is good to go or needs extra work is a slam dunk.

$\Delta G$ kcal/mol	Stability	Fraction denatured
9.6	Small amount of denatured protein	1/10 000 000
8.2		1 / 1 000 000
6.8		1/ 100 000
5.5	Moderate amount of denatured protein	1 / 10 000
4.1		1 / 1 000
2.7	Large amount of denatured protein	1/100
1.3		1/10
0		1/2

## Walk away for hours or days

The only way to get  $\Delta G$  before was a total pain – by hand – and why you probably skipped it altogether. That whole old-school process is completely hands-off with the HUNK. Pick a few conditions or max out on dozens. Add your protein samples, formulations and denaturant. The HUNK does all the work from there and kicks out your data automatically.






## Predict aggregation


Tweak your formulation before long-term stability starts, not after. The HUNK's  $\Delta G_{\text{trend}}$  method gives you the aggregation propensity and pathway. Monitoring  $\Delta G_{\text{trend}}$  over increasing concentrations tells you if your protein will aggregate or not. If it does, you'll know if it happens in the native or denatured state, what fraction is denatured, how much is aggregated, and how much the denatured protein aggregates. Save your SEC time to rubber stamp the winning formulations.

mAb Concentration	$\Delta G$
12 $\mu\text{g/mL}$	8.0
24 $\mu\text{g/mL}$	7.7
48 $\mu\text{g/mL}$	7.5
96 $\mu\text{g/mL}$	7.2
192 $\mu\text{g/mL}$	6.7
384 $\mu\text{g/mL}$	6.5
768 $\mu\text{g/mL}$	6.4



## Add small molecule affinity

When you can't measure small molecule  $K_D$  with your existing platforms – try the HUNK instead.  $\Delta G$ s on the HUNK also let you grab the affinity that small molecule ligands have for receptors. There's no dependence on mass. And because it's a totally solution-based assay, there's no need for immobilization or expensive consumables either.



## Specifications

Description	Specification
Minimum protein concentration per $\Delta G$	25 $\mu\text{g}/\text{mL}$ IgG (protein dependent)
Typical volume required per $\Delta G$	1 mL
Time to measure each $\Delta G$	Average 1 hour
Maximum unattended $\Delta G$ measurements	96
Experiment temperature operating range	Constant 2 °C above ambient to 30 °C (+2 °C)
Protein stock temperature	4 °C to ambient
Fluorescence excitation and emission range	190–890 nm ex, 200–900 nm em photomultiplier
Fluorescence light source	Xe lamp
Fluorescence detection dynamic range	Five orders of magnitude
Environmental conditions	Temperature range: 18–28 °C
Physical	56 cm W x 71 cm D x 87 cm H, 91 kg
Electrical	Auto switching power supply, voltage 100–240 V AC, 50–60 Hz



### Unchained Labs

6940 Koll Center Pkwy, Suite 200  
Pleasanton, CA 94566  
Phone: 1.925.587.9800  
Toll-free: 1.800.815.6384  
Email: [info@unchainedlabs.com](mailto:info@unchainedlabs.com)



Unchained Labs is proudly represented in  
Australia and New Zealand by  
AXT Pty. Ltd.  
1/3 Vuko Pl., Warriewood  
NSW 2102 Australia  
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365  
W. [www.axt.com.au](http://www.axt.com.au) E. [info@axt.com.au](mailto:info@axt.com.au)

© 2016 Unchained Labs. All rights reserved. The  
Unchained Labs logo and HUNK are trademarks  
and/or registered trademarks of Unchained Labs.

Rev B