

LUNATIC



LUNCHAINED
LABS

Get your quant on

Lunatic makes batch quantification of protein, DNA and RNA a no-brainer. All you need is 2 μ L and 5 minutes to measure up to 96 samples. Run them straight-up, even at high concentrations, without ever diluting. Lunatic gets biologics and genomics UV/Vis quantification on the money every time. Just drop, load and read.



96 samples in **5** minutes

Plate format




16 samples in **2** minutes

Chip format



No B.S. workflow

Skip the sample prep, cleaning and worrying about cross-contamination or evaporation – your samples will sit tight for up to 2 hours. Our proprietary Luna Chip consumable has 16 individual microfluidic circuits so you can run up to 16 samples at a time. Or knock out up to 96 on a Luna Plate. So just load and get kick-ass results.





Max out your biologic

Lunatic is the only system out there that can measure biologics at high-throughput and high concentration. It's got dynamic range that covers from 0.02 mg/mL to 150 mg/mL (mAb), so run any protein without ever having to dilute again. Stop the madness of running just one protein at a time – and cleaning up afterwards.




Dig up your genomic's dirt

DNA and RNA samples can be freaking messy. Lunatic analysis software lets you see annoying impurities that classic A260 measurements alone mistake as your DNA or RNA. Find out exactly how much RNA, turbidity, beads and other influencers are screwing up your samples. Ditch the time-sucking, dye-based prep and know right away if your samples are good to go.



Get with the program

GLP labs don't sweat it. Lunatic's software hooks labs up with 21 CFR Part 11 compliant features. We're talking password protection, electronic signatures, full audit trail – the whole package. Lunatic also talks to robots, so you can plug it right in to your current automation workflow.



The screenshot displays the Lunatic Client interface for 21CFR11 compliance. The top navigation bar is green and contains the Lunatic Labs logo, the text "Lunatic Client", the version "21CFR11", and a user profile icon labeled "service".

The "Settings" section includes the following configuration options:

- Force change password after (days): 90
- Maximum inactivity time (minutes): 30
- Maximum number of failed logins: 5
- Minimum characters in password: 5

The "Password needs to contain:" section has the following checked options:

- numbers
- upper case letters
- lower case letters
- special characters

A "SAVE CHANGES" button is located at the bottom right of the settings panel.

The "Audit Trail" section features a table with columns for ID, Date, User, Category, Action, and Details. The table is filtered by "All" for both Category and Action. Two buttons, "EXPERIMENTS INTEGRITY" and "SAVE TO PDF", are positioned above the table.

ID	Date	User	Category	Action	Details
58	2017-06-28 12:01:47	service	User	User login	User login successful (service)
57	2017-06-28 12:01:43		System	Software start	Software start: cDrop
56	2017-06-28 12:01:01		System	Software shutdown	Software shutdown: Lunatic
55	2017-06-28 11:52:35	service	System	Add license	Dongle license added (activation code: "AD284581B072COF5CD751D0E5")
54	2017-06-28 10:24:46	service	User	User login	User login successful (service)
53	2017-06-28 10:24:39		System	Software start	Software start: Lunatic
52	2017-06-27 18:33:02		System	Software shutdown	Software shutdown: cDrop
51	2017-06-27 17:49:43		System	Software shutdown	Software shutdown: Lunatic
50	2017-06-27 17:37:40	service	User	User login	User login successful (service)
49	2017-06-27 17:37:36		System	Software start	Software start: Lunatic
48	2017-06-27 16:04:35	service	Data	Report created	User "service" created a report for experiment done on Mon, Oct 17, 2016
47	2017-06-27 15:54:37	service	Data	New experiment	New experiment done on Mon, Oct 17, 2016 4:19:18 PM, performed on S/
46	2017-06-27 15:25:52	service	User	User login	User login successful (service)

A "BACK" button is located at the bottom left of the interface.

Specifications

Description	Big Lunatic Specifications	Little Lunatic Specifications
Sample size	2 μ L	2 μ L
Number of samples	96 per Luna Plate	16 per Luna Chip
Sample run time	5 minutes per Luna Plate	2 minutes per Luna Chip
IgG concentration range	Luna Plate: 0.02–29 mg/mL High Luna Plate: 0.02–150 mg/mL	Luna Chip: 0.02–29 mg/mL High Luna Chip: 0.02–150 mg/mL
dsDNA concentration range	Luna Plate: 1.5–2,000 ng/ μ L High Luna Plate: 1.5–10,000 ng/ μ L	Luna Chip: 1.5–2,000 ng/ μ L High Luna Chip: 1.5–10,000 ng/ μ L
Path length	Luna Plate: 0.5 mm High Luna Plate: 0.1 and 0.7 mm	Luna Chip: 0.5 mm High Luna Chip: 0.1 and 0.7 mm
Light Source	Xenon flash lamp	Xenon flash lamp
Detector	UV/Vis polychromatic spectrophotometer	UV/Vis polychromatic spectrophotometer
Wavelength range	230–750 nm	230–750 nm
Wavelength accuracy	\leq 0.2 nm	\leq 0.2 nm
Absorbance precision	0.001 OD (0.5 mm path)	0.001 OD (0.5 mm path)
Absorbance accuracy	4% (0.66 OD at 350 nm)	4% (0.66 OD at 350 nm)
Calibration	Integrated self-calibration	Integrated self-calibration
Physical	37 cm W x 46 cm D x 33 cm H, 21 kg	23 cm W x 30 cm D x 28 cm H, 9 kg
Electrical	Universal input, voltage 100–240 V AC, 50–60 Hz	Universal input, voltage 100–240 V AC, 50–60 Hz
Computer	Customer provided computer, Microsoft Windows 7 or later	Preinstalled software, 7" full-color touch screen



Unchained Labs is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au



Unchained Labs

6870 Koll Center Parkway

Pleasanton, CA 94566

Phone: 1.925.587.9800

Toll-free: 1.800.815.6384

Email: info@unchainedlabs.com

© 2017 Unchained Labs. All rights reserved. The Unchained Labs logo, Lunatic and the Lunatic logo are trademarks and/or registered trademarks of Unchained Labs. All other brands or product names mentioned are trademarks owned by their respective organizations.