



SELFRAG ACCESSORIES • SINK STATION

Application

The sink station is used to charge, decant, unload and clean the process vessel.

Mode of operation

The station is equipped with a support frame to hold the process vessel, which can be easily tilted and fixed in a suitable position for decanting, unloading and cleaning. The sink is covered with a perforated metal plate.

The station is equipped with a tap and a handspray attached to a flexible braided hose for easy cleaning of components.

A sedimentation tank with overflow and siphon is located below the sink.



Technical specifications

Design

Galvanized welded steel construction

Outer dimensions

Max. L x W x H	120 x 80 x 140 cm
Working level	91 cm
Weight max.	250 kg
Complete Part-Nr :	202127

Connections:

Water connection	3/4" inch thread
Drain	1" inch thread

