


MODEL 130



Specimen Punch

A precision ground punch and die plate eliminate specimen stress and distortion





MODEL 130

Specimen Punch

The Model 130 Specimen Punch produces high-quality disk specimens from thin metal foils.

- A precision ground punch and die plate eliminate specimen stress and distortion
- A spring-loaded return plunger keeps the disk specimen on the die plate surface for convenient handling
- Available in standard sizes of 1 mm, 2.3 mm, and 3.0 mm
- Other sizes available upon request

DISK SPECIMENS FROM THIN FOILS

The Model 130 Specimen Punch is an ideal accessory to the Fischione Model 110 Automatic Twin-Jet Electropolisher.

The specimen punch is available with a standard base or with a wide base, which provides extra stability when producing thin foils. The specimen punch without a base is shown on the cover, the specimen punch with a base is shown above.



Fischione is proudly represented in Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1.724.325.5444
Fax: +1.724.325.5443
info@fischione.com
www.fischione.com