



MODEL 130

Specimen Punch

The Model 130 Specimen Punch produces high-quality disk specimens from thin metal foils. It is an ideal accessory for the Model 110 Automatic Twin-Jet Electropolisher.

Model 130 Specimen Punch Specifications

Disk size	Available in standard 1 mm, 2.3 mm, and 3 mm sizes Other sizes are available upon request
Base	Available with or without a base
Dimensions	<i>With base:</i> 4.0 in (102 mm) width x 3.1 in (79 mm) height x 4 in (102 mm) depth <i>Without base:</i> 1.0 in (25 mm) width x 2.5 in (64 mm) height x 1.1 in (28 mm) depth
Weight	<i>With base:</i> 0.85 lb (387 g) <i>Without base:</i> 0.25 lb (115 g)
Warranty	One year



Fischione is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1 724.325.5444
Fax: +1 724.325.5443
info@fischione.com
www.fischione.com

©2013 E.A. Fischione Instruments, Inc. All rights reserved.
Document Number SP130 Revision 00 07/2013