



MODEL 160

Specimen Grinder

The Model 160 Specimen Grinder mechanically prethins specimens for transmission electron microscopy (TEM).

Model 160 Specimen Grinder Specifications

Specimen size	Up to 18 mm diam
Dimensions	3.0 in (76 mm) diam x 2.6 in (66 mm) height
Weight	2.4 lb (1.1 kg)
Warranty	One year



Fischione is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au



E.A. Fischione Instruments, Inc.
9003 Corporate Circle
Export, PA 15632 USA
Tel: +1 724.325.5444
Fax: +1 724.325.5443
info@fischione.com
www.fischione.com

©2013 E.A. Fischione Instruments, Inc. All rights reserved.
Document Number SP160 Revision 01 08/2013