



MODEL 9020

Vacuum Pumping Station

Allows for the simultaneous vacuum storage of up to five plasma-cleaned specimens and transmission electron microscope (TEM) specimen holders. Part numbers below are for the Model 9010 Vacuum Storage Container Assembly with Brackets.

Vacuum storage containers may also be purchased separately; contact sales@fischione.com for more information.

Model 9010 Vacuum Storage Container Part Numbers

| | |
|--------------------|---|
| 014-0172-01 | For Zeiss 910, 912, and 922 |
| 014-0172-02 | For JEOL SCSH 1210, 1230, 1400, 2010, 2100, 2200, 2500, 3010, 3200, and ARM200F |
| 014-0173-01 | For FEI Standard: CM12, 20, and 30 |
| 014-0173-02 | For FEI CompuStage CM120, 200, and 300; Tecnai 12, 20, and 30; Titan 60-300 |
| 014-0173-03 | For JEOL SQH 100, 200, 1010, 1011, 1200, and 2000 |
| 014-0173-04 | For JEOL 4000 |
| 014-0173-05 | For Hitachi HF-2000, HD-2000, HD-2300, HD-2700, and H-8100 |
| 014-0173-06 | For Hitachi S-5000, S-5200 and S-5500 |
| 014-0173-07 | Vacuum Storage Container Assembly with Brackets; Hitachi H-9000 |
| 014-0173-08 | For Hitachi HF-3300 |
| Warranty | One year |



E.A. Fischione Instruments, Inc.
 9003 Corporate Circle
 Export, PA 15632 USA
 Tel: +1 724.325.5444
 Fax: +1 724.325.5443
 info@fischione.com
 www.fischione.com



Fischione is proudly represented in Australia and New Zealand by AXT Pty. Ltd.
 1/3 Yuko Pl., Warriewood
 NSW 2102 Australia
 T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
 W. www.axt.com.au E. info@axt.com.au