

## Wizard Classic 1 tubes

1009530

| Tube | Precipitation Reagent                                                  | Buffer                                                                | Salt                      |
|------|------------------------------------------------------------------------|-----------------------------------------------------------------------|---------------------------|
| 1    | 20% (w/v) PEG 8000                                                     | 100 mM CHES/ Sodium hydroxide pH 9.5                                  |                           |
| 2    | 10% (v/v) 2-propanol                                                   | 100 mM HEPES/ Sodium hydroxide pH 7.5                                 | 200 mM Sodium chloride    |
| 3    | 15% (v/v) Reagent alcohol                                              | 100 mM CHES/ Sodium hydroxide pH 9.5                                  |                           |
| 4    | 35% (v/v) MPD                                                          | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Magnesium chloride |
| 5    | 30% (v/v) PEG 400                                                      | 100 mM CAPS/ Sodium hydroxide pH 10.5                                 |                           |
| 6    | 20% (w/v) PEG 3000                                                     | 100 mM Sodium citrate/ Citric acid pH 5.5                             |                           |
| 7    | 10% (w/v) PEG 8000                                                     | 100 mM MES/ Sodium hydroxide pH 6.0                                   | 200 mM Zinc acetate       |
| 8    | 2000 mM Ammonium sulfate                                               | 100 mM Sodium citrate/ Citric acid pH 5.5                             |                           |
| 9    | 1000 mM Ammonium phosphate dibasic                                     | 100 mM Sodium acetate/ Acetic acid pH 4.5                             |                           |
| 10   | 20% (w/v) PEG 2000 MME                                                 | 100 mM Tris base/ Hydrochloric acid pH 7.0                            |                           |
| 11   | 20% (v/v) 1,4-butanediol                                               | 100 mM MES/ Sodium hydroxide pH 6.0                                   | 200 mM Lithium sulfate    |
| 12   | 20% (w/v) PEG 1000                                                     | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Calcium acetate    |
| 13   | 1260 mM Ammonium sulfate                                               | 100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5                    |                           |
| 14   | 1000 mM Sodium citrate tribasic                                        | 100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5                    |                           |
| 15   | 10% (w/v) PEG 3000                                                     | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Lithium sulfate    |
| 16   | 2500 mM Sodium chloride                                                | 100 mM Potassium phosphate monobasic/ Sodium phosphate dibasic pH 6.2 |                           |
| 17   | 30% (w/v) PEG 8000                                                     | 100 mM Sodium acetate/ Acetic acid pH 4.5                             | 200 mM Lithium sulfate    |
| 18   | 1000 mM Potassium sodium tartrate                                      | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Sodium chloride    |
| 19   | 20% (w/v) PEG 1000                                                     | 100 mM Tris base/ Hydrochloric acid pH 7.0                            |                           |
| 20   | 400 mM Sodium phosphate monobasic/ 1600 mM Potassium phosphate dibasic | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Sodium chloride    |
| 21   | 20% (w/v) PEG 8000                                                     | 100 mM HEPES/ Sodium hydroxide pH 7.5                                 |                           |
| 22   | 10% (v/v) 2-propanol                                                   | 100 mM Tris base/ Hydrochloric acid pH 8.5                            |                           |
| 23   | 15% (v/v) Reagent alcohol                                              | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Magnesium chloride |
| 24   | 35% (v/v) MPD                                                          | 100 mM Tris base/ Hydrochloric acid pH 7.0                            | 200 mM Sodium chloride    |
| 25   | 30% (v/v) PEG 400                                                      | 100 mM Tris base/ Hydrochloric acid pH 8.5                            | 200 mM Magnesium chloride |
| 26   | 10% (w/v) PEG 3000                                                     | 100 mM CHES/ Sodium hydroxide pH 9.5                                  |                           |
| 27   | 1200 mM Sodium phosphate monobasic/ 800 mM Potassium phosphate dibasic | 100 mM CAPS/ Sodium hydroxide pH 10.5                                 | 200 mM Lithium sulfate    |
| 28   | 20% (w/v) PEG 3000                                                     | 100 mM HEPES/ Sodium hydroxide pH 7.5                                 | 200 mM Sodium chloride    |
| 29   | 10% (w/v) PEG 8000                                                     | 100 mM CHES/ Sodium hydroxide pH 9.5                                  | 200 mM Sodium chloride    |
| 30   | 1260 mM Ammonium sulfate                                               | 100 mM Sodium acetate/ Acetic acid pH 4.5                             | 200 mM Sodium chloride    |
| 31   | 20% (w/v) PEG 8000                                                     | 100 mM Sodium phosphate dibasic/ Citric acid pH 4.2                   | 200 mM Sodium chloride    |
| 32   | 10% (w/v) PEG 3000                                                     | 100 mM Potassium phosphate monobasic/ Sodium phosphate dibasic pH 6.2 |                           |
| 33   | 2000 mM Ammonium sulfate                                               | 100 mM CAPS/ Sodium hydroxide pH 10.5                                 | 200 mM Lithium sulfate    |
| 34   | 1000 mM Ammonium phosphate dibasic                                     | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            |                           |
| 35   | 20% (v/v) 1,4-butanediol                                               | 100 mM Sodium acetate/ Acetic acid pH 4.5                             |                           |
| 36   | 1000 mM Sodium citrate tribasic                                        | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            |                           |
| 37   | 2500 mM Sodium chloride                                                | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            |                           |
| 38   | 1000 mM Potassium sodium tartrate                                      | 100 mM CHES/ Sodium hydroxide pH 9.5                                  | 200 mM Lithium sulfate    |
| 39   | 20% (w/v) PEG 1000                                                     | 100 mM Sodium phosphate dibasic/ Citric acid pH 4.2                   | 200 mM Lithium sulfate    |
| 40   | 10% (v/v) 2-propanol                                                   | 100 mM MES/ Sodium hydroxide pH 6.0                                   | 200 mM Calcium acetate    |
| 41   | 30% (w/v) PEG 3000                                                     | 100 mM CHES/ Sodium hydroxide pH 9.5                                  |                           |
| 42   | 15% (v/v) Reagent alcohol                                              | 100 mM Tris base/ Hydrochloric acid pH 7.0                            |                           |
| 43   | 35% (v/v) MPD                                                          | 100 mM Potassium phosphate monobasic/ Sodium phosphate dibasic pH 6.2 |                           |
| 44   | 30% (v/v) PEG 400                                                      | 100 mM Sodium acetate/ Acetic acid pH 4.5                             | 200 mM Calcium acetate    |
| 45   | 20% (w/v) PEG 3000                                                     | 100 mM Sodium acetate/ Acetic acid pH 4.5                             |                           |
| 46   | 10% (w/v) PEG 8000                                                     | 100 mM Imidazole/ Hydrochloric acid pH 8.0                            | 200 mM Calcium acetate    |
| 47   | 1260 mM Ammonium sulfate                                               | 100 mM Tris base/ Hydrochloric acid pH 8.5                            | 200 mM Lithium sulfate    |
| 48   | 20% (w/v) PEG 1000                                                     | 100 mM Sodium acetate/ Acetic acid pH 4.5                             | 200 mM Zinc acetate       |



Rigaku is proudly represented in  
Australia and New Zealand by  
AXT Pty. Ltd.  
1/3 Vuko Pl., Warriewood  
NSW 2102 Australia  
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365  
W. [www.axt.com.au](http://www.axt.com.au) E. [info@axt.com.au](mailto:info@axt.com.au)