

Wizard PEG Ion 48 salt screen with PEG 1K and 4K

(96 formulations, 1.7 ml each in a 96-well block plate)

1008655

Well	Precipitation Reagent	Salt (0.2M)
A1	30% (w/v) PEG 1000	Ammonium acetate
A2	30% (w/v) PEG 1000	Ammonium chloride
A3	30% (w/v) PEG 1000	Ammonium fluoride
A4	30% (w/v) PEG 1000	Ammonium formate
A5	30% (w/v) PEG 1000	Ammonium iodide
A6	30% (w/v) PEG 1000	Ammonium nitrate
A7	30% (w/v) PEG 1000	Ammonium phosphate (monobasic)
A8	30% (w/v) PEG 1000	Ammonium phosphate (dibasic)
A9	30% (w/v) PEG 1000	Ammonium sulfate
A10	30% (w/v) PEG 1000	Ammonium sulfite
A11	30% (w/v) PEG 1000	Calcium acetate
A12	30% (w/v) PEG 1000	Calcium chloride
B1	30% (w/v) PEG 1000	Diammonium tartrate
B2	30% (w/v) PEG 1000	Formic acid potassium salt
B3	30% (w/v) PEG 1000	Lithium acetate
B4	30% (w/v) PEG 1000	Lithium chloride
B5	30% (w/v) PEG 1000	Lithium citrate
B6	30% (w/v) PEG 1000	Lithium nitrate
B7	30% (w/v) PEG 1000	Lithium sulfate
B8	30% (w/v) PEG 1000	Magnesium acetate
B9	30% (w/v) PEG 1000	Magnesium chloride
B10	30% (w/v) PEG 1000	Magnesium formate
B11	30% (w/v) PEG 1000	Magnesium nitrate
B12	30% (w/v) PEG 1000	Magnesium sulfate
C1	30% (w/v) PEG 1000	Potassium acetate
C2	30% (w/v) PEG 1000	Potassium chloride
C3	30% (w/v) PEG 1000	Potassium fluoride
C4	30% (w/v) PEG 1000	Potassium iodide
C5	30% (w/v) PEG 1000	Potassium nitrate
C6	30% (w/v) PEG 1000	50mM Potassium phosphate (dibasic)
C7	30% (w/v) PEG 1000	Potassium phosphate (monobasic)
C8	30% (w/v) PEG 1000	Potassium sulfate
C9	30% (w/v) PEG 1000	Potassium thiocyanate
C10	30% (w/v) PEG 1000	Potassium/sodium tartrate
C11	30% (w/v) PEG 1000	Sodium acetate
C12	30% (w/v) PEG 1000	150mM Sodium chloride
D1	30% (w/v) PEG 1000	Sodium citrate
D2	30% (w/v) PEG 1000	Sodium fluoride
D3	30% (w/v) PEG 1000	Sodium formate
D4	30% (w/v) PEG 1000	Sodium iodide
D5	30% (w/v) PEG 1000	Sodium isothiocyanate
D6	30% (w/v) PEG 1000	Sodium nitrate
D7	30% (w/v) PEG 1000	Sodium phosphate (dibasic)
D8	30% (w/v) PEG 1000	Sodium phosphate (monobasic)
D9	30% (w/v) PEG 1000	Sodium sulfate
D10	30% (w/v) PEG 1000	Sodium tartrate (dibasic)
D11	30% (w/v) PEG 1000	Tripotassium citrate
D12	30% (w/v) PEG 1000	Zinc acetate



Rigaku is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au

Well	Precipitation Reagent	Salt (0.2M)
E1	20% (w/v) PEG 4000	Ammonium acetate
E2	20% (w/v) PEG 4000	Ammonium chloride
E3	20% (w/v) PEG 4000	Ammonium fluoride
E4	20% (w/v) PEG 4000	Ammonium formate
E5	20% (w/v) PEG 4000	Ammonium iodide
E6	20% (w/v) PEG 4000	Ammonium nitrate
E7	20% (w/v) PEG 4000	Ammonium phosphate (monobasic)
E8	20% (w/v) PEG 4000	Ammonium phosphate (dibasic)
E9	20% (w/v) PEG 4000	Ammonium sulfate
E10	20% (w/v) PEG 4000	Ammonium sulfite
E11	20% (w/v) PEG 4000	Calcium acetate
E12	20% (w/v) PEG 4000	50mM Calcium chloride
F1	20% (w/v) PEG 4000	Ammonium tartrate dibasic
F2	20% (w/v) PEG 4000	Formic acid potassium salt
F3	20% (w/v) PEG 4000	Lithium acetate
F4	20% (w/v) PEG 4000	Lithium chloride
F5	20% (w/v) PEG 4000	Lithium citrate
F6	20% (w/v) PEG 4000	Lithium nitrate
F7	20% (w/v) PEG 4000	Lithium sulfate
F8	20% (w/v) PEG 4000	Magnesium acetate
F9	20% (w/v) PEG 4000	50mM Magnesium chloride
F10	20% (w/v) PEG 4000	Magnesium formate
F11	20% (w/v) PEG 4000	Magnesium nitrate
F12	20% (w/v) PEG 4000	Magnesium sulfate
G1	20% (w/v) PEG 4000	Potassium acetate
G2	20% (w/v) PEG 4000	100mM Potassium chloride
G3	20% (w/v) PEG 4000	Potassium fluoride
G4	20% (w/v) PEG 4000	Potassium iodide
G5	20% (w/v) PEG 4000	Potassium nitrate
G6	20% (w/v) PEG 4000	Potassium phosphate (dibasic)
G7	20% (w/v) PEG 4000	100mM Potassium phosphate (monobasic)
G8	20% (w/v) PEG 4000	100mM Potassium sulfate
G9	20% (w/v) PEG 4000	Potassium thiocyanate
G10	20% (w/v) PEG 4000	Potassium/Sodium tartrate
G11	20% (w/v) PEG 4000	100mM Sodium acetate
G12	20% (w/v) PEG 4000	Sodium chloride
H1	20% (w/v) PEG 4000	50mM Sodium citrate
H2	20% (w/v) PEG 4000	Sodium fluoride
H3	20% (w/v) PEG 4000	Sodium formate
H4	20% (w/v) PEG 4000	Sodium iodide
H5	20% (w/v) PEG 4000	Sodium isothiocyanate
H6	20% (w/v) PEG 4000	Sodium nitrate
H7	20% (w/v) PEG 4000	Sodium phosphate (dibasic)
H8	20% (w/v) PEG 4000	50mM Sodium phosphate (monobasic)
H9	20% (w/v) PEG 4000	Sodium sulfate
H10	20% (w/v) PEG 4000	Sodium tartrate (dibasic)
H11	20% (w/v) PEG 4000	100mM Potassium citrate tribasic
H12	20% (w/v) PEG 4000	Zinc acetate



Rigaku is proudly represented in
 Australia and New Zealand by
 AXT Pty. Ltd.
 1/3 Vuko Pl., Warriewood
 NSW 2102 Australia
 T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
 W. www.axt.com.au E. info@axt.com.au