

Technical Specifications

	PDB-1	PDB-4	PDB-5-36	PDB-6
dimensions	ca. 150x200 mm	ca. 150x500 mm	ca. 150x750 mm	ca. 15x750 mm
weight	ca. 5 KG	ca. 20 KG	ca. 25 KG	ca. 25 KG
voltage	80-250 V AC 50-60HZ	80-250 V AC 50-60HZ	80-250 V AC 50-60HZ	80-250 V AC 50-60HZ
maximum power	100 W	300 W	600 W	480 W
average power	40 W	150 W	300 W	250 W
temperature range	-10 °C to +100 °C	-10 °C to +100 °C	-10 °C to +100 °C	-10 °C to +100 °C
accuracy	0.1 °C	0.1 °C	0.1 °C	0.1 °C
temperature zones	1	4	5	6
positions per zone	18	18	36	18
hole diameter	10.35 mm	10.35 mm	10.35 mm	10.35 mm
hole depth	80 mm	80 mm	80 mm	80 mm
display	2 color digital	multicolor graphical	multicolor graphical	multicolor graphical
operation	4 keys	touch display	touch display	touch display

All instruments are fulfilling the requirements of the official
AOCS method CD16b-93 Revised Edition 1999



LAIX Technologies is proudly represented in
Australia and New Zealand by
AXT Pty. Ltd.
1/3 Vuko Pl., Warriewood
NSW 2102 Australia
T. +61 (0)2 9450 1359 F. +61 (0)2 9450 1365
W. www.axt.com.au E. info@axt.com.au