



PDB 5-36

The new model in our PDB series: the instruments in our PDB series are the most precise, quiet and energy-efficient. They are operated using a colored touchscreen that provides clearly arranged status information and self-explanatory guidance.

The dry-bathes are available with 1, 4, 5 or 6 temperature zones and can be upgraded with an automatic calibration and an inert gas cover to avoid condensation at low temperatures. All models are ready for use in our automation systems.

Specifications

- -10 °C to +100°C
- 36 positions per zone
- Approx. 50 watt per zone
- 80 to 250 volt AC
- color touch display
- Complies with official AOCS requirements for IUPSC Cd 16b-93

Contact

LAIX Technologies UG
Am Fasanenhang 5
52379 Langerwehe

Tel.: +49 (0) 2409 48798 07
Fax: +49 (0) 2409 48798 08

Email: info@laix-tech.de
Web: www.laix-tech.de



**NEW: now 36 positions
per zone**